

**Amendments to the Specification**

Please replace the paragraph beginning at page 3, line 22, with the following redlined paragraph:

Figure 6 shows nuclear translocation of NF- $\kappa$ B in THP-1 cells (monocyte cell line) untreated (from left, first panel, images; second panel, quantitation of first panel images) and treated with LPS (third panel, images; fourth panel, quantitation of third panel images). Images are from darkfield, NF- $\kappa$ B labeled, brightfield, and 7-AAD nuclear-label labeled.

Please replace the paragraph beginning at page 4, line 6, with the following redlined paragraph:

Figure 11 shows images of nuclear translocation of NF- $\kappa$ B in adherent A-549 cells untreated (from left, first panel, images; second panel, quantitation of first panel images) and treated with IL-1 $\beta$ /TNF- $\alpha$  (third panel, images; fourth panel, quantitation of third panel images). Images are from darkfield, NF- $\kappa$ B labeled, brightfield, and 7-AAD nuclear-label labeled.

Please replace the paragraph beginning at page 14, line 16, with the following redlined paragraph:

By way of background and wishing to be bound by theory, NF- $\kappa$ B resides predominantly in the cytoplasm in resting cells. Activating treatments (*e.g.*, IL-1  $\beta$  /TNF- $\alpha$  or LPS) induce NF- $\kappa$ B translocation into the nucleus in responsive cell types. Thus, the ratio of nuclear to cytoplasmic ~~NF $\kappa$ B~~ NF- $\kappa$ B increases with LPS treatment. Similar to the A-549 cells, NF- $\kappa$ B is translocated from the cytoplasm to the nucleus when the non-adherent human monocyte cell line, THP-1, is exposed to lipopolysaccharide (LPS). Using the identical probing protocol and CCF, again a quantifiable difference in the nuclear localization NF- $\kappa$ B is demonstrated when comparing untreated and LPS-

treated cells (*see* Figures 6 and 9). A nuclear and NF- $\kappa$ B pixel signal correlation analysis CCF was used to quantitate the difference between untranslocated NF- $\kappa$ B and NF- $\kappa$ B translocated to the cell nucleus. The CCF distinguished location-specific (nuclear and cytoplasmic) quantitation of NF- $\kappa$ B to distinguish LPS-treated from untreated THP-1 cells. Thus, the methods of the present disclosure may also be used with non-adherent cells and cell lines.

Please replace the section beginning at page 17, line 1, with the following redlined section:

A. Materials

01. anti-NF $\kappa$ B (F6) : Santa Cruz Biotechnology (Cat. No.SC-8008),  
200  $\mu$ g/ml
02. Alexa Fluor488 donkey anti-mouse IgG: Molecular Probes (~~Cat#~~),  
1.1 mg/ml
03. Streptavidin Alexa Fluor 488: Molecular Probes
04. Recombinant human TNF- $\alpha$  : BD (Cat# 554618, Lot#  
0000056653)
05. Recombinant human IL-1 $\beta$ : ~~ebioscience~~ eBioscience (Cat# 14-  
8018-62, Lot#)
06. A549 cells (ATCC No. CCL-185)
07. Dulbecco's MEM
08. Fetal Calf Serum
09. F-25 Culture Flask
10. 0.25 % trypsin / EDTA
11. Phosphate buffered saline without  $\text{Ca}^{2+}$ / $\text{Mg}^{2+}$  (PBS)
12. 4% PFA/PBS (Fixation Buffer)
13. 0.1% triton X-100/PBS (Perm Buffer)

B. Cell preparation

We used A549 cells cultured in Dulbecco's MEM supplemented with 10% fetal calf serum in an incubator containing 5% CO<sub>2</sub> at 37. A-549 cells were stimulated with or without TNF- $\alpha$  and IL-1 $\beta$  for 45 min to induce nuclear translocation of NF- $\kappa$ B.

01. Culture A549 cells in the T-75 cm<sup>2</sup> culture flask containing 20 ml of the 10% FCS/ Dulbecco's MEM.

02. Stimulate the exponentially growing cells with TNF- $\alpha$  (2.0 ng/ml) and IL-1 $\beta$  (10 pg/ml) for 45 min at 37°C under 5% CO<sub>2</sub> humidified atmosphere.

03. After stimulation, discard media and wash cells with 5-10 ml of PBS.

04. Add 2 ml of 0.25 % trypsin / EDTA to cells, and incubate 37°C for 1 min or until cells have detached.

05. Suspend cells by adding 8 ml of complete DMEM.

06. ~~transfer~~ Transfer the cell suspension to 15 ml centrifuge tube.

07. Centrifuge at 300 x g ~~40'~~ 10', 4°C, and remove media.

08. Fix cells by resuspending @at 10<sup>7</sup> cells/ml in 4% PFA/PBS 30', 4°C.

09. Wash with PBS, then perm cells by resuspending @at 2 x 10<sup>7</sup> cells/ml in 0.1% triton X-100/0.02% EDTA/PBS (Perm) 30', 4°C.

10. Add equal volume of anti-NF $\kappa$ B 20  $\mu$ g/mL in Perm (final mAb concentration of 10  $\mu$ g/mL) 15', 4°C.

11. Wash Perm Buffer.

12. Resuspend 10<sup>7</sup> cells/ml in Perm + AF488 donkey anti-mouse IgG (10  $\mu$ g/mL) 15', 4°C.

13. Filter 70  $\mu$ m mesh and wash with Perm.

14. Resuspend 5 x 10<sup>7</sup> cells/ml Perm + 10  $\mu$ M 7-AAD 5' and run directly on the ImageStream.

## EXAMPLE 2

### INDUCTION OF TRANSLOCATION IN NON-ADHERENT CELLS

Human monocyte cell line THP-1, obtained from ATCC (Rockville, MD), were maintained in RPMI 1640 (Gibco, Grand Island, NY) containing 5% fetal bovine serum, 1 mM sodium pyruvate (Mediatech, Herndon, VA), 100  $\mu$ M nonessential amino acids, 100 U/ml penicillin, 100  $\mu$ g/ml streptomycin, and 2 mM L-glutamine (BioWhittaker, Walkersville, MD) in 5% CO<sub>2</sub> atmosphere at 37°C. The density of exponentially growing cells was less than  $3 \times 10^5$  cells per ml at the time of all treatments. To induce NF- $\kappa$ B translocation into the nucleus from the cytoplasm, cells were treated for 1 hr with LPS.

The following is the experimental procedure for LPS-induced Nuclear Translocation of NF- $\kappa$ B in THP-1 cells.

#### Samples:

1) Unstained and single fluorescent color control samples – start with  $3.0 \times 10^6$  total cells each. In this experiment, controls are:

unstained

NF $\kappa$ B Alexa Fluor488

7-AAD

At the end, resuspend in 100  $\mu$ l 0.1% triton X-100/PBS.

Unstained and NF $\kappa$ B can be mixed and run as one file, then a separate .rif of unlabeled cells can be created in IDEAS. The 7-AAD control must be run separately, because 7-AAD comes off of labeled cells and stains unlabeled cells, confounding compensation. Furthermore, we run the sample with 7-AAD in the buffer to increase staining intensity (washing it away reduces the intensity about four-fold).

2) Experimental samples – start with  $10^7$  total cells for untreated LPS-treated. Stain according to following protocol.

C.A. Materials

- ~~141.~~ anti-NF $\kappa$ B (F6) : Santa Cruz Biotechnology (Cat. No.SC-8008),  
200  $\mu$ g/ml
- ~~152.~~ Alexa Fluor488 donkey anti-mouse IgG: Molecular Probes (~~Cat~~),  
1.1 mg/ml
- ~~163.~~ Streptavidin Alexa Fluor 488: Molecular Probes
- ~~174.~~ Lipopolysaccharide (LPS) from E. Coli 0111B4 : Sigma (Cat#  
L2630, Lot#)
- ~~185.~~ THP-1 cells
- ~~196.~~ RPMI
- ~~207.~~ Fetal Calf Serum
- ~~218.~~ T-75 cm<sup>2</sup> Culture Flask
- ~~229.~~ EDTA
- ~~2310.~~ Phosphate buffered saline without Ca<sup>2+</sup>/Mg<sup>2+</sup> (PBS)
- ~~2411.~~ 4% PFA/PBS (Fixation Buffer)
- ~~2512.~~ 0.1% triton X-100/PBS (Perm Buffer)

D.B. Cell preparation

We used THP-1 cells cultured in RPMI supplemented with 10% fetal calf serum in an incubator containing 5% CO<sub>2</sub> at 37. THP-1 cells were stimulated with or without LPS and for 60 min to induce nuclear translocation of NF- $\kappa$ B.

~~151.~~ Culture THP-1 cells in the T-75 cm<sup>2</sup> culture flask containing 50 ml of the 10% FCS/ RPMI (3x10<sup>5</sup> cells/mL).

~~162.~~ Stimulate the exponentially growing cells with LPS for 60 min at 37°C under 5% CO<sub>2</sub> humidified atmosphere.

~~173.~~ Centrifuge at 300 x g ~~10'~~ 10', 4°C, and remove media.

~~184.~~ Fix cells by resuspending ~~@~~ at 10<sup>7</sup> cells/ml in 4% PFA/PBS 30',  
4°C.

~~195.~~ Wash with PBS, then perm cells by resuspending ~~@~~at  $2 \times 10^7$  cells/ml in 0.1% triton X-100/0.02% EDTA/PBS (Perm) 30', 4°C.

~~206.~~ Add equal volume of anti-NF $\kappa$ B 20  $\mu$ g/mL in Perm (final mAb concentration of 10  $\mu$ g/mL) 15', 4°C.

~~217.~~ Wash Perm Buffer.

~~228.~~ Resuspend  $10^7$  cells/ml in Perm + AF488 donkey anti-mouse IgG (10  $\mu$ g/mL) 15', 4°C.

~~239.~~ Filter 70  $\mu$ m mesh and wash with Perm.

~~2410.~~ Resuspend  $5 \times 10^7$  cells/ml Perm + 10  $\mu$ M 7-AAD 5' and run directly on ImageStream.